

PROTOCOL FOR **GISH** GENOMIC *IN SITU* HYBRIDIZATION

Peng Zhang, Ralf Kynast, and Bernd Friebe
Wheat Genetics Resource Center



Plant Pathology Department
Kansas State University
Manhattan KS 66506-5502
www.ksu.edu/wgrc



TABLE OF CONTENTS

<i>Stock solution recipes</i>	2
<i>Extracting genomic DNA for GISH (CTAB method)</i>	6
<i>Sizing the genomic competitor DNA for GISH</i>	7
<i>F-labeling of genomic DNA for GISH</i>	8
<i>Preparing and pretreating slides for GISH</i>	11
<i>Genomic in situ hybridization (GISH) with F-DNA / competitor DNA</i>	13
<i>Counterstaining (blue or red)</i>	16

STOCK SOLUTION RECIPES.

For all stock solutions to be prepared, use chemicals with highest grade available. All stock solutions must be prepared with deionized distilled H₂O (ddH₂O). For most steps in the DNA handling, it is necessary to autoclave the ddH₂O for 20 min on a liquid cycle in order to destroy any DNase activity and assure sterility. Note, the water quality in these procedures has a stronger effect on the signal-to-noise ratio than the purity grade of the chemicals.

A. Stock solutions stored at room temperature (RT).

- 50% DS** Dissolve 5 g of dextran sulfate in 5 ml ddH₂O.
- 4 M NaOH** Dissolve 16 g of NaOH in ddH₂O to a final volume of 100 ml.
- 1x PBS pH 7.4** Dilute 26 ml of 5 M NaCl + 20 ml of 500 mM PB pH 7.4 with 954 ml ddH₂O to a final volume of 1,000 ml.
- 2 x SSC pH 7.0** Dilute 100 ml of 20 x SSC pH 7.0 with 900 ml ddH₂O to a final volume of 1,000 ml.
- 1x TAE buffer pH 8.0** Dilute 20 ml of 50 x TAE buffer pH 8.0 with 980 ml ddH₂O to a final volume of 1,000 ml.
- 1.5x extraction buffer (1L) (for CTAB (hexadecyltrimethylammonium bromide) method)**
75 ml 1M Tris, pH 8.0 + 15 g CTAB + 150 ml 0.5M EDTA, pH 8.0 + 0.29715 g 1,10 phenanthroline + 61.3575 g NaCl + (add before use) 1.5 ml 2-mercaptoethanol. Add ddH₂O to a final volume of 1,000 ml.

B. Stock solutions stored at 4°C (refrigerator).

- 1% acetocarmine** Dissolve 10 g of carmine (Fisher) in 1,000 ml of 45% glacial acetic acid. Add boil-ezers and reflux (in hood) for 24 hr. Filter into a dark bottle.
- Antifade** Use Vectashield (Vector Laboratories, H-1000) for mounting slides.
- 100 mM citric acid** Dissolve 21.0 g of citric acid in ddH₂O to 1,000 ml.
- FA** Use formamide with high grade (e.g., Merck or Fisher).
- 50% FA in 2x SSC** Dilute 50 ml of formamide + 10 ml of 20 x SSC with 40 ml ddH₂O.
- LB-amp broth, pH 7.0 (1% Bacto tryptone, 0.5% Bacto yeast extract, 1% NaCl, 50 µg/ml ampicillin), sterile**
Dissolve 4 g of Bacto-tryptone (Difco) + 2 g of Bacto-yeast extract (Difco) + 4 g of NaCl in 380 ml ddH₂O. Titrate pH to 7.0 with NaOH and bring volume to 400 ml with ddH₂O. Autoclave for 20 min at 121°C, cool to room temperature, add 400 µl of 50 mg/ml ampicillin under sterile conditions.
- Mcllvaine buffer pH 7.0** Dilute 17.65 ml of 100 mM citric acid + 32.94 ml of 500 mM Na₂HPO₄ with 49.41 ml ddH₂O
- 1 M MgCl₂** Dissolve 20.33 g of MgCl₂·6H₂O in ddH₂O to a final volume of 100 ml.
- 3 M Na Ac pH 5.2** Dissolve 40.81 g of sodium acetate·3H₂O in 30 ml ddH₂O, titrate pH to 5.2 with glacial acetic acid, and dilute with ddH₂O to a final volume of 100 ml.
- 5 M NaCl** Dissolve 29.22 g of NaCl in ddH₂O to a final volume of 100 ml.

GENOMIC IN SITU HYBRIDIZATION

- 500 mM NaH₂PO₄** Dissolve 68.99 g of NaH₂PO₄·H₂O in ddH₂O to a final volume of 1,000 ml.
- 500 mM Na₂HPO₄** Dissolve 88.99 g of Na₂HPO₄·2H₂O in ddH₂O to a final volume of 1,000 ml.
- 500 mM PB (phosphate buffer) pH 7.4 (350 mM Na₂HPO₄, 150 mM NaH₂PO₄)**
Add 500 mM Na₂HPO₄ + 500 mM NaH₂PO₄ to pH 7.4.
- 10x PBS (phosphate buffered saline) pH 7.4 (1.3 M NaCl, 70 mM Na₂HPO₄, 30mM NaH₂PO₄)**
Dilute 26 ml of 5 M NaCl + 20 ml of 500 mM PB pH 7.4 with 54 ml ddH₂O.
- PFA** Use high-grade paraformaldehyde (e.g., Merck or Fisher).
- 20 x SSC (saline sodium citrate) pH 7.0 (3 M NaCl, 0.3 M Na₃ citrate)**
Dissolve 175.32 g of NaCl + 88.23 g of Na₃ citrate·2H₂O in 850 ml ddH₂O, titrate pH to 7.0 with NaOH/HCl, add ddH₂O to a final volume of 1,000 ml.
- 50 x TAE buffer pH 8.0 (2 M Tris acetate pH 8.0, 50 mM Na₂EDTA)**
Dissolve 18.6 g of Na₂EDTA + 242.24 g of Tris base in 700 ml ddH₂O, titrate pH to 8.0 with glacial acetic acid (~ 57 ml), add ddH₂O to a final volume of 1,000 ml.
- 500 mM Na₂EDTA pH 8.0** Dissolve 18.61 g of disodium ethylene diamine tetraacetate·2H₂O + 2 g of NaOH in 80 ml ddH₂O by vigorous stirring. Titrate pH to 8.0 with NaOH.
- 1 M Tris.Cl buffer pH 7.5** Dissolve 121.12 g of Tris base in 800 ml ddH₂O, titrate pH to 7.5 with concentrated HCl (~ 65 ml), and add ddH₂O to a final volume of 1,000 ml.
- 1 M Tris.Cl buffer pH 8.0** Dissolve 121.12 g of Tris base in 800 ml ddH₂O, titrate pH to 8.0 with concentrated HCl (~ 42 ml), and add ddH₂O to a final volume of 1,000 ml.
- TE buffer pH 8.0 (10 mM Tris.Cl pH 8.0, 1 mM Na₂EDTA)**
Dilute 20 µl of 1 M Tris-Cl pH 8.0 + 4 µl of 500 mM Na₂EDTA pH 8.0 with 1,976 µl ddH₂O.
- 15 mM NH₄Ac in 80% EtOH (for CTAB method)**
Take 1 ml of 7.5M stock of NH₄Ac. Use 80% EtOH to a final volume of 500 ml.

C. Stock solutions stored at -20°C (freezer).

- 10 µg/µl carrier DNA** Use sonicated, denatured DNA from salmon testes (Sigma D-7656) and store in aliquots.

Counterstains

- 100 µg/ml DAPI [Excitation_{max}: 344 nm, Emission_{max}: 449 + 488 nm]**
Dissolve 1 mg of 4',6-diamidino-2-phenylindole·2HCl in 10 ml ddH₂O and store in 1-ml aliquots.
- 100 µg/ml PI [Excitation_{max}: 340 + 530 nm, Emission_{max}: 615 nm]**
Dissolve 1 mg of propidium iodide in 10 ml ddH₂O and store in 1-ml aliquots.
- 5 U/µl DNA Pol I** Use *E. coli* DNA Polymerase I (Kornberg enzyme) (Boehringer # 642711) and store in 10-µl aliquots.
- 30 U/µl DNase I** Use pancreatic DNase I solution (e.g. Boehringer # 776785) and store in 10-µl aliquots.

G E N O M I C I N S I T U H Y B R I D I Z A T I O N

10x NT buffer (Nick Translation buffer) pH 7.5 (500 mM Tris.Cl pH 7.5, 50 mM MgCl₂, 0.5 mg/ml BSA)

Dilute 1 ml of 1 M Tris.Cl pH 7.5 + 100 µl of 1 M MgCl₂ + 10 µl of 10% BSA with 890 µl ddH₂O and store in 100-µl aliquots.

100 mM DTT

Dissolve 155 mg of 1,4-dithiothreitol + 33 µl of 3 M Na-acetate pH 5.2 in ddH₂O to 10 ml and store in 1-ml aliquots.

1.3 mM dNTP mix for nick translation (400 µM dATP, 400 µM dGTP, 400 µM dCTP, 100 µM dTTP)

Use 100 mM solutions of the lithium salts of the dNTPs (Boehringer). Dilute 4 µl of dATP + 4 µl of dGTP + 4 µl of dCTP + 1 µl of dTTP with 987 µl ddH₂O and store in 100-µl aliquots.

1 mM F-x-dUTP

Use 1 mM solution of an appropriate Fluorophore-x-dUTP (Amersham, Boehringer, DuPont, Enzo, or Molecular Probes; see Table on p. 5).

1% RNase A in 10 mM Tris.Cl pH 7.5, 15 mM NaCl (DNase free)

Dissolve 10 mg of RNase A (Sigma) + 10 µl of 1 M Tris.Cl pH 7.5 + 3 µl of NaCl in 987 µl ddH₂O, incubate in boiling waterbath for 15 min, cool slowly, and store in aliquots.

Alternative counterstains

Counterstain	Catalog #	Excitation _{max}		Emission _{max}	
		nm	color	nm	color
POPO-1	Molecular Probes P-3580	434	Violet	456	Blue
DAPI	Molecular Probes D-1306	359	Ultraviolet	461	Blue
Hoechst 33342	Molecular Probes H-3570	355	Ultraviolet	465	Blue
Hoechst 33258	Molecular Probes H-3569	365	Ultraviolet	480	Blue
BOBO-1	Molecular Probes B-3582	462	Blue	481	Blue
YOYO-1	Molecular Probes Y-3601	491	Blue	509	Green
SYTOX Green	Molecular Probes S-7020	504	Green	523	Green
TOTO-1	Molecular Probes T-3600	514	Green	533	Green
Lucifer Yellow VS	Sigma L-3510	428	Violet	540	Green
SYTO 15	Molecular Probes S-7577	515	Green	546	Yellow
Acridine Yellow	Sigma A-2520	470	Blue	550	Yellow
POPO-3	Molecular Probes P-3584	534	Green	570	Orange
BOBO-3	Molecular Probes B-3586	570	Orange	604	Red
Etbr	Molecular Probes E-3565	518	Green	605	Red
PI	Molecular Probes P-1304	535	Green	617	Dark red
YOYO-3	Molecular Probes Y-3606	612	Dark red	631	Dark red
TOTO-3	Molecular Probes T-3604	642	Dark red	660	Far red

Alternative F-x-dUTPs

F-x-dUTP	Catalog #	Excitation _{max}		Emission _{max}	
		nm	color	nm	color
ChromaTide Cascade Blue-7-dUTP	Molecular Probes C-7612	400	Violet	420	Blue
FluoroBlue Coumarin-4-dUTP	Amersham RPN 2123	350	Ultraviolet	440	Blue
Coumarin-5-dUTP	DuPont NEL-415	402	Violet	443	Blue
Amino-methylcoumarin-6-dUTP	Boehringer 1 534 386	342	Ultraviolet	450	Bright blue
Amino-methylcoumarin-6-dUTP	Enzo 42851	342	Ultraviolet	450	Bright blue
Amino-methylcoumarin-13-dUTP	Enzo 42852	342	Ultraviolet	450	Bright blue
ChromaTide BODIPY FL-14-dUTP	Molecular Probes C-7614	505	Green	513	Green
Fluorescein-12-dUTP	Enzo 42831	495	Blue	517	Green
Fluorescein-12-dUTP	DuPont NEL-413	494	Blue	517	Green
Fluorescein-12-dUTP	Boehringer 1 373 242	494	Blue	517	Green
ChromaTide Fluorescein-12-dUTP	Molecular Probes C-7604	494	Blue	518	Green
ChromaTide Alexa 488-5-dUTP	Molecular Probes C-11397	490	Blue	519	Green
FluoroGreen Fluorescein-11-dUTP	Amersham RPN 2121	490	Blue	520	Green
ChromaTide Oregon Green 488-5-dUTP	Molecular Probes C-7630	496	Blue	524	Green
ChromaTide Rhodamine Green-5-dUTP	Molecular Probes C-7629	502	Green	527	Green
ChromaTide Alexa 532-5-dUTP	Molecular Probes C-11398	525	Green	550	Yellow
Tetramethylrhodamine-6-dUTP	DuPont NEL-416	550	Yellow	570	Orange
ChromaTide Tetramethylrhodamine-6-dUTP	Molecular Probes C-7606	550	Yellow	570	Orange
ChromaTide Alexa 546-5-dUTP	Molecular Probes C-11401	555	Yellow	570	Orange
FluoroLink Cy3-dUTP	Amersham PA 53022	550	Yellow	570	Orange
ChromaTide BODIPY TMR-14-dUTP	Molecular Probes C-7616	542	Green	574	Orange
FluoroRed Rhodamine-4-dUTP	Amersham RPN 2122	545	Yellow	575	Red
Tetramethyl-rhodamine-6-dUTP	Boehringer 1 534 378	550	Yellow	575	Red
Tetramethylrhodamine-5-dUTP	Enzo 42841	551	Yellow	575	Red
ChromaTide Tetramethylrhodamine-5-dUTP	Molecular Probes C-7606	555	Yellow	580	Red
Lissamine-5-dUTP	DuPont NEL-418	570	Orange	588	Red
ChromaTide Alexa 568-5-dUTP	Molecular Probes C-11399	575	Red	600	Red
Texas Red-dUTP	DuPont NEL-417	593	Red	612	Dark red
ChromaTide Alexa 594-5-dUTP	Molecular Probes C-11400	590	Red	615	Dark red
ChromaTide Texas Red-12-dUTP	Molecular Probes C-7631	595	Red	615	Dark red
ChromaTide Texas Red-5-dUTP	Molecular Probes C-7608	595	Red	615	Dark red
ChromaTide BODIPY TR-14-dUTP	Molecular Probes C-7618	589	Red	617	Dark red
ChromaTide BODIBY 630/650-14-dUTP	Molecular Probes C-11395	632	Dark red	648	Dark red
ChromaTide BODIPY 650/665-14-dUTP	Molecular Probes C-11396	650	Dark red	668	Far red
FluoroLink Cy5-dUTP	Amersham PA 55022	649	Dark red	670	Far red

EXTRACTING GENOMIC DNA FOR GISH (CTAB METHOD)**A. Isolating DNA.**

1. Grind 2–5 g of green leaves to a fine powder in a mortar and pestle with liquid nitrogen.
2. Transfer the powder to a 50-ml tube containing 20 ml of 1.5x extraction buffer. Mix thoroughly but gently with a precooled spatula. Put in a water bath at 65°C for 3 hr; gently mix every 30 min.
3. Cool tube to RT and add chloroform:isoamylalcohol (24:1) to the total volume of 40 ml. Mix thoroughly but gently by hand for 30 min until the top layer is milky green in color.
4. Centrifuge at 4,000 rpm for 10 min at RT.
5. Transfer the supernatant with a disposable plastic pipette to a new tube and add 2 ml (1/10 volume) of 3M NaAc and an equal volume (the volume of the suspension) of isopropanol (–20°C). Mix by hand for 10 min. Leave at RT for 10 min.
6. Using a 1-ml tip to retain the DNA clot, remove the isopropanol solution. Wash the clot with 15 mM NH₄Ac in 80% EtOH (the total volume is 35 ml). Leave at RT for 10 min, mixing occasionally.
7. Drain the ethanol solution and dry for 10–30 min. Add 650 µl TE depending upon the size of the DNA pellet. Then add 1 µl RNase (10 mg/ml). Leave at 4°C overnight. If DNA is needed immediately, it may be incubated at 37°C for 1 hr.

B. Purifying the DNA.

1. Make sure that the DNA is completely dissolved. Add 1:1 phenol:chloroform equal to the total volume of TE added in step 7 (250 µl each for 500 µl TE). Mix thoroughly but gently for 5–10 min and centrifuge at 13,500 rpm for 5–7 min at RT. Larger amounts of DNA will require a longer time.
2. Transfer the upper phase to a new tube. Extract with an equal volume of chloroform:isoamyl alcohol (24:1) (500 µl). Mix well and centrifuge for 1 min.
3. Transfer the supernatant to a new tube. Extract with an equal volume of chloroform:isoamylalcohol (24:1) again. Mix well and centrifuge for 1 min.
4. Add 1/10 volume (65 µl) of 3M NaAc, pH 5.2 and 1 volume of isopropanol (650 µl). Mix thoroughly but gently. Leave at RT for 10 min. Use a glass hook to keep the DNA clot.
5. Drain the solution and fill tube twice with cold 70% ethanol (–20°C) and once with RT 70% ethanol to clean the DNA. Leave at RT for 10 min.
6. Drain ethanol and dry the pellet as much as possible. Resuspend pellet in an amount of TE depending upon the DNA yield (usually 500 µl).

C. Minigeelectrophoresis.

1. Mix the samples (5 µl):
 - a. 1 µl of genomic DNA + 3 µl ddH₂O + 1 µl of 5x DNA loading buffer and
 - b. 5 µl of DNA-ladder with DNA-loading buffer.
2. Load samples into the wells and electrophoresis.
3. Check gel on a UV light transilluminator for the band of genomic DNA.

SIZING GENOMIC COMPETITOR DNA FOR GISH

Genomic DNA to be used as the unlabeled competitor DNA is sized to an average length of ~ 200 nucleotide pairs in order to enable the GISH process to occur in highly defined stringency conditions and to promote the best penetration into the condensed structure of a mitotic or meiotic plant chromosome or chromatin of an interphase nucleus.

A. Pretreating the DNA.

1. Dilute the extracted genomic DNA with TE buffer pH 8.0 to a concentration ~1 µg/µl
2. Autoclave the diluted DNA for 10 min at 121°C, 15 PSI.
3. Cool the autoclaved DNA and store at -20°C in 100-µl aliquots.

B. Minigeelectrophoresis.

1. Mix the samples (5 µl):
 - a. 1 µl of 1/10 diluted untreated genomic DNA (100 ng) + 3 µl ddH₂O + 1 µl of 5x DNA loading buffer and
 - b. 5 µl of 1-kb DNA ladder with DNA-loading buffer.
2. Load samples into the wells and electrophoresis.
3. Check gel on a UV light transilluminator for the band of the competitor DNA at about 200 bp.

F-LABELING OF GENOMIC DNA FOR GISH

Genomic DNA is labeled by nick translation in the presence of F-x-dUTP. Nick translation (Rigby *et al.* 1977: J. Mol. Biol. 113, 237-251) is based on the introduction of random single-stranded nicks in double-stranded DNA by pancreatic DNase I. The three activities of the *E. coli* DNA polymerase I then catalyze the addition of nucleotide residues to the 3'-hydroxyl terminus of a nick (5'→3' polymerase activity) with the simultaneous elimination of nucleotides from the 5'-phosphoryl terminus (5'→3' exonuclease activity). As nucleotides are removed and new ones added, the nick is linearly moved, or translated along the DNA strand. In the presence of F-x-dUTP, a dTTP analog, about 50 or 60 % of the dTTP residues are replaced with F-x-dUTP. Using a low reaction temperature (15°C) avoids synthesis of snap back DNA structures by optimal proofreading of the DNA polymerase I (3'→5' exonuclease activity). After denaturation, the final probe length (ssDNA) should be ~ 400 nucleotides for best penetration of the labeled probe into the condensed structure of a mitotic or meiotic plant chromosome or chromatin of an interphase nucleus.

A. Pretreating the DNA.

1. Shear DNA for 200x to a medium length of \pm 12 or 15 kbp with a sterile 1-ml syringe.
2. Dilute the sheared DNA with ddH₂O to a concentration of ~1 μ g/ μ l.

B. Diluting DNase I.

1. **Immediately prior to use**, remove an aliquot of DNase I stock from -20°C freezer and put on ice.
2. Dilute (1/500*) the thawed DNase I aliquot on ice: 1 μ l of DNase I (30 U/ μ l) + 499 μ l ddH₂O

* The dilution must effect a nicking activity for a final probe length of 400 nucleotides (ssDNA), i.e., about 600 or 700 basepairs (dsDNA). Depending on the purity and medium length of the cloned DNA and the enzyme activity of the DNase I stock solution, the needed dilution can vary from about 1/50 (for BACs or YACs) to about 1/2,000 (for microcloned DNA of about 2.5 kbp).

→ 1/50 dilution	»	600 mU/ μ l	→ 1/500 dilution	»	60 mU/ μ l
→ 1/100 dilution	»	300 mU/ μ l	→ 1/1,000 dilution	»	30 mU/ μ l
→ 1/300 dilution	»	100 mU/ μ l	→ 1/2,000 dilution	»	15 mU/ μ l

C. Nick-translation.

1. Set up the reaction mix (50 μ l) in a sterile 200- μ l ml safe-lock tube on ice.

Component	Volume (μ l)	Mass per reaction	Final concentration
ddH ₂ O	30-x		
10x NT buffer	5	2.5 μ mol Tris.Cl pH 7.5	50 mM
		250 nmol MgCl ₂	5 mM
		2.5 μ g BSA	50 ng/ μ l
100 mM DTT	5	500 nmol	10 mM
1.3 mM dNTP mix, equimolar ^(*)	5	2 nmol dATP	40 μ M
		2 nmol dGTP	40 μ M
		2 nmol dCTP	40 μ M
		500 pmol dTTP	10 μ M
1 mM F-x-dUTP, equimolar *	2	2 nmol	40 μ M
DNA	x	1 μ g	20 ng/ μ l
DNase I, diluted	1	60 mU	1.2 mU/ μ l
DNA Pol I	2	10 U	0.2 U/ μ l

* 2.5 nmol (dTTP + F-x-dUTP) for mol fraction T takes into account the > 25% lower incorporation rate of fluorochrome labeled dUTPs into the DNA than that of unmodified dTTP.

2. Thoroughly, but carefully, mix the solutions and incubate for 2 h at 15°C.
3. Stop incubation by placing the tube on ice, then incubate at 65°C for 7 min to inactivate the enzymes.
4. Purify the probe using a QIAquick Nucleotide Removal Kit (QIAGEN 28304).
 - a. Add 500 μ l Buffer PN to the reaction mixture and mix.
 - b. Place a QIAquick spin column in a 2-ml collection tube.
 - c. To bind DNA, apply the sample to the QIAquick column and centrifuge for 1 min at 6,000 rpm.
 - d. Discard the flow-through and place QIAquick column back into the same tube.
 - e. Wash QIAquick column by adding 750 μ l of Buffer PE and centrifuge for 1 min at 6,000 rpm.
 - f. Discard the flow-through and place the QIAquick column back in the same tube, which should be empty. Centrifuge for an additional 1 min at 10,000 x g.
 - g. Place the QIAquick column in a clean 1.5-ml microcentrifuge tube.
 - h. Elute DNA by adding 50 μ l of Buffer EB to the center of the QIAquick membrane, let the column stand for 1 min, and then centrifuge for 1 min at 10,000 x g.
5. Store the probe at -20°C

D. Minigelectrophoresis.

1. Mix samples (5 μ l):
 - a. 2 μ l of F-probe (~40 ng) + 2 μ l ddH₂O + 1 μ l of 5x DNA loading buffer and
 - b. 5 μ l of DNA 1-kb ladder with loading buffer.
2. Load samples into the wells and electrophorese in 1x TAE buffer pH 8.0 for 1 h at RT.
3. Check the gel on a UV light transilluminator for the band of the probe at about 600 to 700 bp.

E. Spot test.

1. Spot 1 μ l of F-probe (~20 ng) onto a small sheet (2 cm x 3 cm) of nylon membrane and air dry for ~10 min.
2. Examine fluorescence intensity under light at appropriate wavelength.

PREPARING AND PRETREATING SLIDES FOR GISH

Among plant tissues containing actively dividing cells, root-tip meristems are the most commonly used for making mitotic chromosome preparations. Other tissues may be used, such as leaf meristems, calli, or protoplasts. Root-tip meristem cells may be obtained from germinating seeds, newly formed roots in soil or hydroponics, or root tips from tissue-culture-grown plants. Root tip collection includes a pretreatment in order to arrest as many metaphase cells as possible. The common method is the so-called ice-water treatment followed by the fixation in modified Carnoy's fixative (absolute ethanol : glacial acetic acid = 3:1).

A. Preparing the slides.

1. Germinate seeds in distilled water on filter paper in a petridish at RT for 2–3 days (2 cm long roots).
2. Cut roots and pretreat in distilled water in tubes on ice for 24 h.
3. Fix roots in fresh fixative (3 volumes of 100% EtOH + 1 volume of glacial acetic acid) at RT.
4. Store roots at 4°C until ready for use.
5. Put roots in 1% acetocarmine for a few min or until the tip of the root is dark red.
6. Cut the very tip of the root and gently squeeze out meristem cells from the root tip onto a slide.
7. Squash meristem cells in a drop of 45% acetic acid under a cover slip (18 x 18 mm). Examine mitotic index under the phase-contrast microscope.
8. Freeze the slide upside down on dry ice. Remove cover slip with the edge of a double-sided razor blade and incubate in 100% EtOH for 5 min at RT.
9. Air dry.

B. RNase A treatment.

1. **Immediately prior to use**, remove an 8- μ l aliquot of RNase A stock from the –20°C freezer and thaw on ice.
2. Dilute (1/100) the thawed RNase A aliquot on ice: 8 μ l of 1% RNase A + 80 μ l of 20 x SSC + 712 μ l ddH₂O.
3. Add 100 μ l of diluted RNase A onto one slide, cover with a plastic cover slip (20 x 20 mm), and incubate in a wet chamber for 45 min at 37°C.
4. Wash slides 3x by gently shaking in 2x SSC for 5 min each at RT.

C. Postfixation.

1. **Immediately prior to use**, prepare 100 ml of 4 % PFA:
 - Warm 80 ml ddH₂O to 60–70°C
 - Add 4 g of PFA
 - Add 100 µl of 4 M NaOH
 - Add 10 ml of 10 x PBS
 - Cool to RT and add ddH₂O to 100 ml.
2. Incubate slides in PFA for 10 min at RT in a fume hood.
3. Wash slides 3x in 2 x SSC for 5 min each at RT.

D. Dehydration.

1. Incubate slides in a 70%, 95%, and 100% EtOH series, for 3 min each at RT.
2. Air dry slides at RT.

GENOMIC IN SITU HYBRIDIZATION (GISH) WITH F-DNA/COMPETITOR DNA

A. Denaturing slides.

1. Set the temperature of a hot plate to 80°C.
2. Denature slides in 100 µl of 70% FA / 2 x SSC for 2 min at 80°C.
3. Dehydrate slides in 70% EtOH (-20°C) for 5 min, 95% EtOH (RT) for 3 min, and 100% EtOH (RT) for 3 min.
4. Air dry slides and prepare for GISH.

B. Fluorescence in situ hybridization.

1. Prepare 30 µl of hybridization solution (HS) per slide, which contains:

5 µg of autoclaved salmon testes DNA as carrier DNA (~ 150 bp*)
 X µg of autoclaved genomic DNA as unlabeled competitor DNA (~ 200 bp*)
 40–60 ng of nick translated genomic DNA as fluorochrome labeled probe DNA (~ 400 bp*)

* The amount of competitor DNA depends on the distance of the relationship between the plant species used as the probe and competitor DNAs (See table below).

DNA amount competitor(s) : probe	Competitor plant species (unlabeled DNA)	Probe plant species (F-labeled DNA)
25 : 1	<i>Aegilops umbellulata</i> (U)	<i>Ae. comosa</i> subsp. <i>comosa</i> (M)
25 : 1	<i>Ae. comosa</i> subsp. <i>comosa</i> (M)	<i>Ae. umbellulata</i> (U)
30 : 1	<i>Triticum aestivum</i> (ABD)	<i>S. cereale</i> (R ^{cer})
40 : 1	<i>T. aestivum</i> (ABD)	<i>Ae. markgrafii</i> (C)
45 : 1	<i>T. aestivum</i> (ABD)	<i>Ae. kotschyi</i> (US)
45 : 1	<i>T. aestivum</i> (ABD)	<i>Ae. ventricosa</i> (DN)
45 : 1	<i>T. aestivum</i> (ABD)	<i>Ae. comosa</i> subsp. <i>comosa</i> (M)
45 : 1	<i>T. aestivum</i> (ABD)	<i>Ae. umbellulata</i> (U)
45 : 1	<i>Hordeum vulgare</i> (H ^{vul})	<i>H. bulbosum</i> (H ^{bul})
50 : 1	<i>Secale strictum</i> subsp. <i>strictum</i> (R ^{mon})	<i>S. cereale</i> (R ^{cer})
50 : 1	<i>S. vavilovii</i> (R ^{vav})	<i>S. cereale</i> (R ^{cer})
50 : 1	<i>S. sylvestre</i> (R ^{syl})	<i>S. cereale</i> (R ^{cer})
xx : 1	<i>Ae. speltooides</i> (S)	<i>T. urartu</i> (A)
xx : 1	<i>Ae. tauschii</i> (D)	<i>T. urartu</i> (A)
xx : 1 + yy : 1	<i>Ae. speltooides</i> (S) + <i>Ae. tauschii</i> (D)	<i>T. urartu</i> (A)
xx : 1	<i>T. urartu</i> (A)	<i>Ae. speltooides</i> (S)
xx : 1	<i>Ae. tauschii</i> (D)	<i>Ae. speltooides</i> (S)
xx : 1 + yy : 1	<i>T. urartu</i> (A) + <i>Ae. tauschii</i> (D)	<i>Ae. speltooides</i> (S)
xx : 1	<i>T. urartu</i> (A)	<i>Ae. tauschii</i> (D)
xx : 1	<i>Ae. speltooides</i> (S)	<i>Ae. tauschii</i> (D)
xx : 1 + yy : 1	<i>T. urartu</i> (A) + <i>Ae. speltooides</i> (S)	<i>Ae. tauschii</i> (D)
xx : 1	<i>T. turgidum</i> subsp. <i>durum</i> (AB)	<i>Ae. tauschii</i> (D)
xx : 1	<i>T. urartu</i> (A)	<i>T. turgidum</i> subsp. <i>durum</i> (AB)
xx : 1	<i>Ae. speltooides</i> (S)	<i>T. turgidum</i> subsp. <i>durum</i> (AB)
xx : 1	<i>Ae. tauschii</i> (D)	<i>T. turgidum</i> subsp. <i>durum</i> (AB)
120 : 1	<i>Trifolium repens</i> (NO)	<i>Tr. ambiguum</i> L. (A)
120 : 1	<i>Tr. ambiguum</i> (A)	<i>Tr. repens</i> L. (NO)

G E N O M I C I N S I T U H Y B R I D I Z A T I O N

for 4 slides use 120 μ l of HS in a sterile 1.5-ml save-lock tube

Component	Volume (μ l)	Final concentration
ddH ₂ O	12.2	
dFA	60.0	50 %
20x SSC	9.0	1.5x SSC
50 % DS	24.0	10 %
carrier DNA (10 μ g/ μ l)	2.0	167 ng/ μ l autoclaved salmon testes DNA (~150 bp*)
competitor DNA (1 μ g/ μ l)	4.8	4.8 μ g autoclaved genomic DNA (~200 bp, e.g. wheat*)
probe DNA (20 ng/ μ l)	8.0	160 ng F-labeled genomic DNA (~400 bp, e.g. rye)

* GISH stringencies (%) according the formula **stringency (%) = 100 (%) - f_m(T_m (°C) - T_h (°C))**.

2. Boil the HS for 5 min.
3. Chill on ice for 5 min and spin down condensed water.
4. Add 30 μ l of HS onto each slide and cover with a plastic cover slip (20 x 20 mm).
5. Incubate in a wet chamber overnight at 37°C.

C. Post-hybridization washes.

1. Wash slides 2x by gently shaking in 2 x SSC for 5 min at 42°C.
2. Wash slides 1x by gently shaking in 50% FA in 2 x SSC for 10 min at 42°C.
3. Wash slides 2x by gentle shaking in 2 x SSC for 5 min at 42°C.
4. Wash slides in 2 x SSC for 5 min at RT in the dark.

G E N O M I C I N S I T U H Y B R I D I Z A T I O N

$$T_m = \text{melting temperature } (^{\circ}\text{C}) = 81.5 (^{\circ}\text{C}) + 16.6 \log_{10} M + 0.41 (\% \text{ G+C}) - 500/n - 0.61 (\% \text{ formamide})$$

M: 2x SSC = 0.3900 M Na⁺
 1.5x SSC = 0.2925 M Na⁺
 1x SSC = 0.1950 M Na⁺
 0.5x SSC = 0.0975 M Na⁺

(% G+C) = 45.5 %
 n ~ 400 (probe)
 n ~ 200 (competitor)
 n ~ 150 (carrier)

(% formamide) = 50 %
 (% formamide) = 40 %
 (% formamide) = 30 %
 on the assumption that f_m = mismatch factor = 1

$$\text{stringency } (\%) = 100 + \text{hybridization temperature} - \text{melting temperature} = 100 + T_h - T_m$$

For the following tables: PR = probe DNA, CO = competitor DNA, CA = carrier DNA, and D = denatured.

T _h (°C)	50% FA/2x SSC			50% FA/1.5x SSC			50% FA/1x SSC			50% FA/0.5x SSC		
	PR	CO	CA	PR	CO	CA	PR	CO	CA	PR	CO	CA
37	75 %	77 %	77 %	77 %	79 %	79 %	80 %	82 %	82 %	85 %	87 %	87 %
42	80	82	82	82	84	84	85	87	87	90	92	95
47	85	87	87	87	89	89	90	92	92	95	97	97
52	90	92	92	92	94	94	95	97	97	100	D	D
57	95	97	97	97	99	99	100	D	D	D	D	D
62	100	D	D	D	D	D	D	D	D	D	D	D
67	D	D	D	D	D	D	D	D	D	D	D	D
72	D	D	D	D	D	D	D	D	D	D	D	D

T _h (°C)	40% FA/2x SSC			40% FA/1.5x SSC			40% FA/1x SSC			40% FA/0.5x SSC		
	PR	CO	CA	PR	CO	CA	PR	CO	CA	PR	CO	CA
37	69 %	71 %	71 %	71 %	73 %	73 %	74 %	76 %	76 %	79 %	81 %	81 %
42	74	76	76	76	78	78	79	81	81	84	86	86
47	79	81	81	81	83	83	84	86	86	89	91	91
52	84	86	86	86	88	88	89	91	91	94	96	96
57	89	91	91	91	93	93	94	96	96	99	D	D
62	94	96	96	96	98	98	99	D	D	D	D	D
67	99	D	D	D	D	D	D	D	D	D	D	D
72	D	D	D	D	D	D	D	D	D	D	D	D

T _h (°C)	30% FA/2x SSC			30% FA/1.5x SSC			30% FA/1x SSC			30% FA/0.5x SSC		
	PR	CO	CA	PR	CO	CA	PR	CO	CA	PR	CO	CA
37	63 %	65 %	65 %	65 %	67 %	67 %	68 %	70 %	70 %	73 %	75 %	75 %
42	68	70	70	70	72	72	73	75	75	78	80	80
47	73	75	75	75	77	77	78	80	80	83	85	85
52	78	80	80	80	82	82	83	85	85	88	90	90
57	83	85	85	85	87	87	88	90	90	93	95	95
62	88	90	90	90	92	92	93	95	95	98	100	100
67	93	95	95	95	97	97	98	100	100	D	D	D
72	98	100	100	100	D	D	D	D	D	D	D	D

COUNTERSTAINING (BLUE OR RED)

DAPI: Excitation_{max}: 344 nm (**ultraviolet**) Emission_{max}: 449 nm (**blue**) + 488 nm (**blue**)

PI: Excitation_{max}: 340 nm (**ultraviolet**) + 530 nm (**green**) Emission_{max}: 615 nm (**red**)

for 4 slides: 2 ml of counterstaining solution

DAPI counterstain (blue).

1. **Immediately prior to use**, remove an aliquot of DAPI stock from a -20°C freezer and thaw on ice.
2. Dilute (1/250) DAPI on ice: 8 μl of DAPI (100 ng/ μl) + 1,992 μl of McIlvaine buffer.
3. Add 500 μl of diluted DAPI onto each slide and incubate for 2 min at RT.
4. Wash slides with ~ 3 ml of McIlvaine buffer from a pipette
5. Mount slides with a small drop of Vectashield antifade and cover with cover glass (24 x 30 mm).
6. Store slides in the dark at 4°C .

PI counterstain (red).

1. **Immediately prior to use**, remove an aliquot of PI stock from a -20°C freezer and thaw on ice.
2. Dilute (1/2,500) PI on ice: 0.8 μl of PI (100 ng/ μl) + 2,000 μl of 1 x PBS.
3. Add 500 μl of diluted PI onto each slide and incubate for 2 min at RT.
4. Wash slides with ~ 3 ml of 1 x PBS from a pipette.
5. Mount slides with a small drop of Vectashield antifade and cover with cover glass (24 x 30 mm).
6. Store slides in the dark at 4°C .